

CORE Prosthetic KIT

CORE Prosthetic KIT Specifications

Order Code : CPK

• A Universal Kit which is applicable for all kinds of Abutment Screw normally used in dentistry field

1.2 Hex Driver (Ratchet/Engine)

Used for tightening Abutment Screw



1.25 Hex Driver (Ratchet/Engine)

Used for tightening Abutment Screw



0.9 Hex Driver (Ratchet/Engine)

Used for tightening Abutment Screw




0.5 Slot Driver (Ratchet)

Used for when it is not available to connect with torque due to screw damaged



SCS Driver (Ratchet/Engine)

Used for tightening Abutment Screw



Unigrip Driver (Ratchet/Engine)

Used for tightening Abutment Screw



1.2/1.25 Angled Driver (Ratchet)

Used for tightening Abutment Screw which the change of angle is needed



Bottom Plate

Torque Wrench

Used for connecting with torque

