

CORE Sinus Master KIT Specifications

Order Code : CSMK

- Integrated Crestal & Lateral Approach
- Available for Guided/Non-guided Type

Crestal Guide Drill

Used for drilling



Initial Guide Drill

Used for first perforation



Tissue Punch

Used for removing gingivitis



Lateral Guide Drill

Used for lateral guided serratus



Lateral Drill

Used for lateral lift




Lateral Wall Expander

Used for expanding the lateral wall perforated





Depth Gauge

Used for checking the depth perforated



Membrane Lifter

Used for water pressure augmentation



Non-Guided Stopper

Select and use according to the depth perforated



Guided Stopper

Select and use according to the depth perforated



Bone Pusher

Used for injecting bone in drill hole



Curette 2

Used for serrating and separating Membrane



Curette 1

Used for serrating and separating Membrane



Bone Carrier

Used for injecting bone

